

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

DISINFECTION APPLICATIONS

Hard Non-Porous Surface Disinfection: To Clean, Disinfect and Deodorize Hard, Non-Porous Surfaces: For heavily soiled areas, a preliminary cleaning is required. Apply, wipe, spray, or dip Danolyte® at 500 ppm FAC to hard, non-porous surfaces with a cloth, wipe, mop or sponge. Treated surfaces must remain wet for 10 minutes. Allow surfaces to air dry. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical devices prior to sterilization or high-level disinfection.

Special Instructions for Cleaning Prior to Disinfection against Clostridium difficile endospores:

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks, or eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with clean cloth, mop, and/or sponge saturated with product intended for disinfection. Cleaning should include vigorous wiping and scrubbing, until visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Materials Disposal: Cleaning materials used that may contain feces/wastes should be disposed of immediately in accordance with local regulations for infectious materials disposal.

Killing Clostridium difficile: Clean hard, non-porous surfaces by removing gross filth, loose dirt, debris, blood/bodily fluids, etc. Apply this product and let stand for 10 minutes.

Special Instructions for Using this product to Clean and Decontaminate against HIV on Surfaces/Objects Soiled with Blood/Body Fluids: This product kills HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care settings (e.g. hospitals, nursing homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) associated with AIDS.

KILL THE GERMS,
NOT THE PLANET

danolyte

Hospital Level Disinfectant

disinfectant
+ sanitizer



Personal Protection: When handling items soiled with blood or body fluids, use appropriate barrier protection such as disposable latex gloves, gowns, masks, and eye coverings.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product.

Contact Time: Apply to area to be treated. Let stand for 10 minutes. Cleaning materials used that may contain feces/wastes should be disposed of immediately in accordance with local regulations for infectious materials disposal.

Disposal of Infectious Material: Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

Storage: Store in a closed dark plastic container away from direct sunlight. Store container in a cool dry area. Product may be disposed in a sanitary sewer.

Pesticide Disposal: Wastes resulting from the use of this product must be disposed of on-site or at an approved waste disposal facility.

Container Disposal: Refillable container. Refill this container with same product only. Do not reuse this container for any other purpose. Cleaning before refilling is the responsibility of the refiller. Cleaning the container before final disposal is the responsibility of the person disposing the container. To clean the container before final disposal, empty the remaining contents from this container into application equipment or mix tank. Fill the container about 10 percent full with water. Agitate vigorously or recirculate water with the pump for two minutes. Repeat this rinsing procedure two more times. Then offer for recycling if available, or puncture and dispose of in a sanitary landfill, or by incineration, or by other procedures allowed by state and local authorities.

Do not use on steel, aluminum, silver, or chipped enamel. Prolonged contact with metal may cause pitting or discoloration. First test in an inconspicuous place for color washout or contact incompatibility.

EPA Reg. No. 91582-1

EPA Est. No. 91582-KS-3

danolyte

KILL THE GERMS. NOT THE PLANET.

DATE PRODUCED:

Manufactured by:

Danolyte Global LLC.
10601 W.79th Street
Shawnee, KS 66214



EPA REGISTERED CLEAN SCENT NO ALCOHOL

NO HARSH CHEMICALS
EFFECTIVE SANITIZER FOR
FOOD CONTACT SURFACES

ACTIVE INGREDIENTS:

Hypochlorous Acid _____ 0.046%

OTHER INGREDIENTS: _____ 99.954%

TOTAL _____ 100.000%

CONTAINS 500PPM (FAC)

KEEP OUT OF REACH OF CHILDREN

(See first aid statement on back panel)



Environmental Commitment

This product rapidly breaks down entirely to salt water. Not harmful to septic and sewer systems. This bottle is coded for recyclers. Check to see if recycling facilities accept colored HDPE. Contains no phosphorous, VOCs, alcohol or phenols. Low Odor. Fresh and clean scent.

- ✓ This product was tested according to AOAC efficacy testing requirements
- ✓ Tested according to AOAC germicidal spray method for hospital disinfectant
- ✓ Meets surface disinfection recommendations under OSHA's Bloodborne Pathogen Guidelines or Standards.

FIRST AID

Have the product container or label with you when calling a poison control center or doctor or going for treatment. For general information on product use, etc, call the National Pesticide Information Center (NPIC) at 1-800-858-7378. For emergencies, call the poison control center at 1-800-222-1222.

NET CONTENTS: 1 US Quart/0.95 Liters

